



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand® Dispensette S Organic, Analog, DE-M,
2.5-25 ml, without RDV**

**€448.70
plus VAT &
Shipping**

Product Images



Description

The Dispensette® S bottle-top dispenser is the number 1 for dispensing directly from the supply bottle. It is equipped with everything that makes dispensing liquids safer and pleasantly easy. Thanks to the special operating principle of the piston, the Dispensette® S is smooth and effortless to operate during liquid aspiration and dispensing. The volume can be quickly and easily adjusted and securely fixed. The Dispensette® S bottle-top dispenser is proven in demanding continuous use and in the use of aggressive media. Dispensette® S Organic (color code yellow): Ideal for dispensing organic solvents, e.g. chlorinated and fluorinated hydrocarbons such as trichlorotrifluoroethane and dichloromethane, or acids such as concentrated HCl and HNO₃ (except HF), as well as for trifluoroacetic acid (TFA), tetrahydrofuran (THF) and peroxides. Analog version: Analog volume display, quickly adjustable dosing volume, securely fixed by internal toothed bar. Easy adjustment with adjustment key.

Additional Information

No.	BN-1366
Manufacturer (Brand)	Brand
VGKL number	104943150
Manufacturer EAN	4033378482273
Transport temperature	Room temperature

