



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Brand® Dispensette S Organic, Digital, DE-M,  
1-10 ml, without RDV**

**€396.70**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

The Dispensette® S bottle-top dispenser is the number 1 for dispensing directly from the supply bottle. It is equipped with everything that makes dispensing liquids safer and pleasantly easy. Thanks to the special operating principle of the piston, the Dispensette® S is smooth and effortless to operate during liquid aspiration and dispensing. The volume can be quickly and easily adjusted and securely fixed. The Dispensette® S bottle-top dispenser is proven in demanding continuous use and in the use of aggressive media. Dispensette® S Organic (color code yellow): Ideal for dispensing organic solvents, e.g. chlorinated and fluorinated hydrocarbons such as trichlorotrifluoroethane and dichloromethane, or acids such as concentrated HCl and HNO<sub>3</sub> (except HF), as well as for trifluoroacetic acid (TFA), tetrahydrofuran (THF) and peroxides. Digital version: Particularly easy-to-read digital display. The selected dosing volume can be set precisely and reproducibly. Easy Calibration technology for adjustment in seconds without tools. Automatic marking when the factory setting is changed.

## Additional Information

---

No.	BN-1378
Manufacturer (Brand)	Brand
VGKL number	104943340
Manufacturer EAN	4033378482136
Transport temperature	Room temperature

