



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand® Dispensette S, Analog, DE-M, 0.1-1 ml, with RDV

**€395.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Dispensette® S bottle-top dispenser is the number 1 for dispensing directly from the supply bottle. It is equipped with everything that makes dispensing liquids safer and pleasantly easy. Thanks to the special operating principle of the piston, the Dispensette® S is smooth and effortless to operate during liquid aspiration and dispensing. The volume can be quickly and easily adjusted and securely fixed. The Dispensette® S bottle-top dispenser is proven in demanding continuous use and in the use of aggressive media. Dispensette® S (color code red): Ideal for dispensing aggressive reagents, e.g. concentrated alkalis and acids such as HPO<sub>4</sub>, H<sub>2</sub>SO<sub>4</sub> (except e.g. HCl, HNO<sub>3</sub> and HF), salt solutions and a variety of organic solvents. Analog version: Analog volume display, quickly adjustable dosing volume, securely fixed by internal toothed bar. Easy adjustment with adjustment key.

## Additional Information

---

No.	BN-1321
Manufacturer (Brand)	Brand
VGKL number	104946101
Manufacturer EAN	4033378482402
Transport temperature	Room temperature

