



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand® Dispensette S, Digital, DE-M, 0.1-1 ml, without RDV

**€396.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

The Dispensette® S bottle-top dispenser is the number 1 for dispensing directly from the supply bottle. It is equipped with everything that makes dispensing liquids safer and pleasantly easy. Thanks to the special operating principle of the piston, the Dispensette® S is smooth and effortless to operate during liquid aspiration and dispensing. The volume can be quickly and easily adjusted and securely fixed. The Dispensette® S bottle-top dispenser is proven in demanding continuous use and in the use of aggressive media. Dispensette® S (color code red): Ideal for dispensing aggressive reagents, e.g. concentrated alkalis and acids such as HPO4, H2SO4 (except e.g. HCl, HNO3 and HF), salt solutions and a wide range of organic solvents. Digital version: Particularly easy-to-read digital display. The selected dosing volume can be set precisely and reproducibly. Easy Calibration technology for adjustment in seconds without tools. Automatic marking when the factory setting is changed.

## Additional Information

No.	BN-1342
Manufacturer (Brand)	Brand
VGKL number	104946310
Manufacturer EAN	4033378487612
Transport temperature	Room temperature

