



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hirschmann® opus dispenser , 10 ml, 100-240 Volt, EU, US, UK plug

€2,366.00
plus VAT &
Shipping

Product Images



Description

opus® dispenser, 10 ml, 100-240 Volt (EU, US, UK plug), electronic motorized bottle top dispenser with touch screen module "dispenser", with reagent feedback system, with declaration of conformity according to DIN 12600, serial or USB interface, complete with screw thread A 45 mm, 3 thread adapters (A 32, A 38 and S 40), 1 suction hose with union nut, 1 dosing ejector unit and 1 power supply unit, Data Power cable 0.75 m, operating instructions, 1 valve mounting lock, with individual serial number and calibration certificate.

Additional Information

No.	NH-1623
Manufacturer (Brand)	Hirschmann
VGKL number	103781010
Transport temperature	Room temperature

