



Thermo Scientific F1-ClipTip multichannel pipettes

Product Images



Description

F1-ClipTip multichannel pipettes Rely on ClipTip technology that securely locks into place to ensure a complete seal of each channel with minimal effort to pick up and drop tips. Available in 8- and 12-channel models. Compatible tips: ClipTip 12.5 Ext Description: F1-ClipTip Number of channels: 8 Volume (metric): 1 - 10 micro;l increments: 0.02 micro;l Color Code: Pink Unit Size: Each

Additional Information

No.	PO-0198
Manufacturer (Brand)	Thermo Scientific

