



## Thermo Scientific F1-ClipTip multichannel pipettes

### Product Images

---



## Description

F1-ClipTip multichannel pipettes Rely on ClipTip technology that securely locks into place to ensure a complete seal of each channel with minimal effort to pick up and drop tips. Available in 8- and 12-channel models. Compatible tips: ClipTip 12.5 Ext Description: F1-ClipTip Number of channels: 12 Volume (metric): 1 - 10 micro;l increments: 0.02 micro;l Color Code: Pink Unit Size: Each

## Additional Information

|                      |                   |
|----------------------|-------------------|
| No.                  | PO-0199           |
| Manufacturer (Brand) | Thermo Scientific |

