



Thermo Scientific F1-ClipTip multichannel pipettes

Product Images



Description

F1-ClipTip multichannel pipettes Rely on ClipTip technology that securely locks into place to ensure a complete seal of each channel with minimal effort to pick up and drop tips. Available in 8- and 12-channel models. Compatible tips: ClipTip 200 Description: F1-ClipTip Number of channels: 12 Volume (metric): 10 to 100 micro;l increments: 0.2 micro; l Color Code: Yellow Unit Size: Each

Additional Information

No.	PO-0195
Manufacturer (Brand)	Thermo Scientific

