



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex® Socorex Qualitix Colorless macro tip, 10 mL, autoclavable, 40x100/bag

€1,074.84
plus VAT &
Shipping

Product Images

Description

Socorex Qualitix® Colorless macro tip, 10 mL, autoclavable, 40x100/bag

Additional Information

No.	SR-0557
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	10 ml
Volume MIN	1 ml
sterile	No

