



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex® Socorex Qualitix Colorless macro tip, 2 mL, autoclavable, 30x250/bag

€1,071.72
plus VAT &
Shipping

Product Images

Description

Socorex Qualitix® Colorless macro tip, 2 mL, autoclavable, 30x250/bag

Additional Information

No.	SR-0553
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	2000 µl
Volume MIN	200 µl
sterile	No

