



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Socorex® Socorex Qualitix Colorless tip, 1000 µL, autoclavable, oblique graduated S

**€78.83**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Socorex Qualitix® Colorless Tip, 1000 µL, autoclavable, oblique graduated tip, 6x96/rack

## Additional Information

---

No.	SR-0551
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	No

