




Thermo Scientific ClipTip™ pipette tips

Product Images



Description

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 50, non-sterile Product line: ClipTip Volume (metric): 50 µl Filtered: No Autoclavable: autoclavable Length (imperial): 2.01 in. Length (metric): 5.1 cm Compatibility: Thermo Scientifictrade; F1-ClipTip Compatibility: Thermo Scientifictrade; F1-ClipTip Format: Rack Tip Shape: Low Retention Packaging: 10 x 96 tips/rack Color Code: Purple Unit Size: Case of 960

Additional Information

No.	PO-0419
Manufacturer (Brand)	Thermo Scientific

