



Thermo Scientific ClipTip™ pipette tips

Product Images



Description

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 20, refill, non-sterile Product line: ClipTip Volume (metric): 20 µl; I Filtered: No Autoclavable: Non-autoclavable Length (angle): 1.8 inches Length (metric): 4.6 cm Compatibility: Thermo Scientifictrade; F1-ClipTip Compatibility: Thermo Scientifictrade; F1-ClipTip Format: Reload Tip Shape: Low Retention Packaging: 10 x 96 tips/insert (reload tower) Color Code: Pink Unit Size: Case of 960

Additional Information

No.	PO-0416
Manufacturer (Brand)	Thermo Scientific

