



Thermo Scientific ClipTip™ pipette tips

Product Images



Description

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 1000, refill, non-sterile Product line: ClipTip Volume (metric): 1000 µl Filtered: No Autoclavable: Non-autoclavable Length (metric): 3.74 in. Length (metric): 9.5 cm Compatibility: Thermo Scientifictrade; E1-ClipTip Compatibility: Thermo Scientifictrade; F1-ClipTip Format: Reload Tip Shape: Low Retention Packaging: 8 x 96 tips/insert (reload tower) Color Code: Blue Unit Size: Case of 768

Additional Information

No.	PO-0439
Manufacturer (Brand)	Thermo Scientific

