




## Thermo Scientific ClipTip™ pipette tips

### Product Images

---



### Description

---

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 50, refill, non-sterile Product line: ClipTip Volume (metric): 50 µl; I Filtered: No Autoclavable: Not autoclavable Length (angle): 2.01 in. Length (metric): 5.1 cm Compatibility: Thermo Scientifictrade; F1-ClipTip Compatibility: Thermo Scientifictrade; F1-ClipTip Format: Refill Tip Shape: Low Retention Packaging: 10 x 96 tips/insert Color Code: Purple Unit Size: Case of 960

# Additional Information

|                      |                   |
|----------------------|-------------------|
| No.                  | PO-0421           |
| Manufacturer (Brand) | Thermo Scientific |

