



Thermo Scientific ClipTip™ pipette tips

Product Images



Description

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 300 Ext, reloading tower, sterile Product line: ClipTip Volume (metric): 300 µl Filtered: Without filter Autoclavable: Non-autoclavable Length (anglo): 3.94 in Compatibility: Thermo Scientifictrade; E1-ClipTip, Thermo Scientifictrade; F1-ClipTip Format: refill Tip Shape: low retention Packaging: 10 x 96 tips/insert (reload tower) Color Code: Orange Unit Size: Each

Additional Information

No.	PO-0435
Manufacturer (Brand)	Thermo Scientific

