



Thermo Scientific ClipTip™ pipette tips

Product Images



Description

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 1250, sterile Product line: ClipTip Volume (metric): 1250 µl; I Filtered: No Length (anglo): = 4.06 in. Length (metric): 10.3 cm Compatibility: Thermo Scientifictrade; E1-ClipTip Compatibility: Thermo Scientifictrade; F1-ClipTip Format: Rack, 8 x 96 Tip Shape: Low Retention Packaging: 8 x 96 tips/rack Color Code: Turquoise Unit Size: Case of 768

Additional Information

No.	PO-0442
Manufacturer (Brand)	Thermo Scientific

