



Thermo Scientific ClipTip™ pipette tips

Product Images



Description

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 200, refill, non-sterile Product line: ClipTip Volume (metric): 200 µl; I Filtered: No Autoclavable: Not autoclavable Length (angle): 2.2 in. Length (metric): 5.6 cm Compatibility: Thermo Scientifictrade; E1-ClipTip Compatibility: Thermo Scientifictrade; F1-ClipTip Format: Refill Tip Shape: Extra Long Packaging: 10 x 96 tips/insert Color Code: Yellow Unit Size: Case of 960

Additional Information

No.	PO-0425
Manufacturer (Brand)	Thermo Scientific

