



Thermo Scientific Finntip pipette tips

Product Images

Description

Finntip pipette tips Ensure accurate and precise pipetting. Designed with all our liquid handling experience, these tips have smooth, hydrophobic surfaces that prevent liquid retention. Description: Finntip 250 Universal, 0.5 - 250 µl, 5.2 cm, yellow For use with (devices): Thermo Scientific Finnpipette 250 Product line: Finntip Volume (metric): 250 µl; I Filtered: No Length (angle): 2 inch Length (metric): 5.2 cm Certification/Compliance: Manufactured to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD mark) Compatibility: Thermo Scientifictrade; Finnpipettetrade; 250 Format: Box Tip shape: Universal fit Packaging: 500 tips Unit Size: Case of 500

Additional Information

No.	PO-0264
Manufacturer (Brand)	Thermo Scientific

