



Thermo Scientific Labtip pipette tips

Product Images

Description

Labtip pipette tips Standard pipette tips for general liquid transfer requirements. Description: Labtip 200, Hinged For use with (devices): Thermo Scientifictrade; Finnpipettetrade; Product line: labtip Volume (metric): 200 µL Filtered: Non-filtered Compatibility: Thermo Scientifictrade; Finnpipettetrade; Format: Rack Unit Size: Case of 3,840

Additional Information

No.	PO-0397
Manufacturer (Brand)	Thermo Scientific

