



Nalgene® CryoBox™ freezer racks

Product Images

Description

CryoBox™ Freezer Racks Conveniently store standard CryoBoxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Vertical Dimensions (L x W x H): 5.625 x 5.5 x 16 in. (14.3 x 14 x 40.6 cm) Capacity: 4 CryoBoxes/Rack Type: Freezer Rack Unit Size: Each

Additional Information

No.	CD-2165
Manufacturer (Brand)	Nalgene
VGKL number	238750372

