



Nalgene® CryoBox™ freezer racks

Product Images

Description

CryoBox™ Freezer Racks Conveniently store standard cryoboxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Vertical Dimensions (L x W x H): 3.375 x 3.25 x 19.75 in. (8.4 x 8.3 x 50.2 cm) Capacity: 9 CryoBoxes/Rack Type: Freezer Rack Unit Size: Each

Additional Information

No.	CD-2164
Manufacturer (Brand)	Nalgene
VGKL number	238750352

