



Nalgene® CryoBox™ freezer racks

Product Images

Description

CryoBox™ Freezer Racks Conveniently store standard cryoboxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Vertical Dimensions (L x W x H): 5.625 x 5.5 x 27.75 in. (14.3 x 14 x 70.5 cm) Capacity: 7 CryoBoxes/Rack Material: Stainless Steel Type: Freezer Rack Unit Size: Each

Additional Information

No.	CD-2166
Manufacturer (Brand)	Nalgene
VGKL number	238750373

