



# Nalgene® CryoBox™ freezer racks

## Product Images

## Description

CryoBox™ Freezer Racks Conveniently store standard cryoboxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Vertical Dimensions (L x W x H): 5.625 x 5.5 x 8.75 in. (14.3. x 14 x 22.5 cm) Capacity: 4 CryoBoxes/Rack Quantity: 4/cs. Type: Cryogenic Box Storage Rack, Vertical Unit Size: Case of 4

## Additional Information

|                      |           |
|----------------------|-----------|
| No.                  | CD-1537   |
| Manufacturer (Brand) | Nalgene   |
| VGKL number          | 238750361 |

