



Nalgene® CryoBox™ freezer racks

Product Images

Description

CryoBox™ Freezer Racks Conveniently store standard cryoboxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Vertical Dimensions (L x W x H): 5.625 x 5.5 x 19.75 in. (14.3 x 14 x 50.2 cm) Capacity: 9 CryoBoxes/Rack Quantity: 4/cs. Type: Cryogenic Box Storage Rack, Vertical Unit Size: Case of 4

Additional Information

No.	CD-1538
Manufacturer (Brand)	Nalgene
VGKL number	238750362

