



Nalgene® CryoBox™ freezer racks

Product Images

Description

CryoBox™ Freezer Racks Conveniently store standard cryoboxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Horizontal Dimensions (L x W x H): 5.625 x 8.875 x 23.25 in. (22.5 x 59.0 x 14.3 cm) Capacity: 16 CryoBoxes/Rack Quantity: Each Type: Cryogenic Box Storage Rack, Horizontal Unit Size: Case of 2

Additional Information

No.	CD-1540
Manufacturer (Brand)	Nalgene
VGKL number	238750382

