



Nalgene® polycarbonate test tube racks

Product Images



Description

Nalgene Polycarbonate Test Tube Racks Thermo Scientific Nalgene polycarbonate test tube racks feature handles for easy transport.
Description: Test tube racks Array: 4 x 10 Autoclavable Capacity: Tubes (13 to 16 mm) Positions: 40 Shape: Rectangle Stackable: Yes
Type: Test Tube Rack For use with (devices): 13 -16mm tubes Number per box: 4 Number per package: 1 Length (Anglo-American): 11 inches Width (anglo): 4.5 inches Diameter (metric) Opening: 16 mm Unit Size: Case of 4

Additional Information

No.	CD-1616
Manufacturer (Brand)	Nalgene

