



Nalgene® polycarbonate test tube racks

Product Images



Description

Nalgene Polycarbonate Test Tube Racks Thermo Scientific Nalgene polycarbonate test tube racks feature handles for easy transport. Description: Test tube racks Array: 4 x 10 Autoclavable Capacity: Tubes (16 to 20 mm) Positions: 40 Shape: Rectangle Stackable: Yes Type: Test Tube Rack For use with (devices): 16 -20mm tubes Number per box: 4 Number per package: 1 Length (Anglo-American): 11 inches Width (anglo): 4.5 inches Diameter (metric) Opening: 20 mm Unit Size: Case of 4

Additional Information

No.	CD-1617
Manufacturer (Brand)	Nalgene

