



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® polypropylene test tube and draining racks

€38.68 As low
as
€31.24
plus VAT &
Shipping

Product Images



Description

Nalgene Polypropylene Test Tube and Draining Rack Thermo Scientific Nalgene polypropylene test tube and draining racks are ideal for test tubes, reaction tubes as well as test tubes. Description: Tube hook rack Array: 8 x 12 Autoclavable: Autoclavable Capacity: Tubes (10 to 13 mm) Positions: 96 Shape: Rectangle Type: test tube plug-in rack For use with (devices): Tubes 10 to 13 mm Number per box: 8 Number per package: 2 Dimensions (L x W x H): 17.9 x 12.7 x 6.4 cm Length (metric): 17.9 cm Length (anglo): 7 in. Width (metric): 12.7 cm Width (Anglo): 5 in. Height (metric): 6.4 cm Height (Anglo): 2.5 in. Diameter (metric) Opening: 13 mm Unit Size: Case of 8

Additional Information

No.	CD-1690
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

