



Nalgene® Pipette Vials

Product Images

Description

Nalgene pipette jars You can insert or sterilize glass pipettes in these jars. The lid indicates evaporation of the cleaning solution and protects the contents from contamination, while the bottom flange increases stability. For use with: Pipettes up to 20.3 cm (8 in.) high (Anglo-American): 20.25 in. Height (metric): 51.4 cm Size: A Description: Pipet Jar XLPE D Diameter (Anglo-American): 6.5 in. Diameter (metric): 16.5 cm Storage Conditions: Store the storage jar only in a container, sink, or draining surface large enough to hold the contents of the jar in the event of a spill. Quantity per package: 1 Quantity per case: 6 Unit Size: Case of 6

Additional Information

No.	CD-1561
Manufacturer (Brand)	Nalgene
VGKL number	248752424

