



## Thermo Scientific Finnpiquette™ Stand

### Product Images

---



## Description

Finnpipette™ Stands Thermo Scientific Finnpi­pette pipette stands provide safe and convenient storage for manual and electronic pipettes. Description: F-Stand, Straight Stand, White Capacity: Manual pipettes (6) Positions: 6 Type: Pipette Stand Unit Size: Each

## Additional Information

|                      |                   |
|----------------------|-------------------|
| No.                  | PO-0392           |
| Manufacturer (Brand) | Thermo Scientific |

