



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY® eye emergency station BR326495 in dust-tight stainless steel cabinet V4A

**€219.65
plus VAT &
Shipping**

Product Images



Description

Including pictogram on the door. With one pH Neutral emergency rinsing bottle à 1000 ml phosphate buffer solution (4.9%) and one eye rinsing bottle à 1000 ml sterile sodium chloride solution (0.9%). Particularly suitable for use in workplaces where chemicals and also foreign bodies can get into the eye and a very long rinsing time is required. - Dimensions wall cabinet: (H x W x D): 355 x 265 x 100 mm - Durability: 3 years - Rinsing time: approx. 20 minutes - According to DIN EN 15154-4:2009

Additional Information

No.	FG-0055
Manufacturer (Brand)	B-SAFETY
Transport temperature	Room temperature

