



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY® eye wash station BR324005 with 2 eye wash bottles

€112.37
plus VAT &
Shipping

Product Images



Description

in a dust-protected wall cabinet with viewing window, material sheet steel, white powder-coated. Including pictogram on the door. With two eye wash bottles of 500 ml sterile sodium chloride solution (0.9%). Particularly suitable for use in workplaces where dust and dirt accumulate (e.g. workshops). Very easy and quick to open in case of first aid. - Dimensions wall cabinet (H x W x D): 355 x 265 x 100 mm - Shelf life: 3 years - Rinsing time: approx. 10 minutes - According to DIN EN 15154-4:2009

Additional Information

No.	FG-0047
Manufacturer (Brand)	B-SAFETY
Transport temperature	Room temperature

