



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY® Plum iBox 2 with 1 x 200 ml pH Neutral and 1 x 500 ml Eye Wash

€57.63
plus VAT &
Shipping

Product Images



Description

Eye emergency station in dust-protected wall box with 1 x 200 ml pH Neutral and 1 x 500 ml eye rinsing solution and separate pictogram panel with mirror. Particularly suitable for use in workplaces where acids and alkalis as well as foreign bodies can enter the eye and where dust and dirt accumulate. Product features - dimensions (H x W x D): 425 x 294 x 135 mm - weight: 1.26 kg - the eye wash bottles are clearly visible - easy and clear recognition of the expiration date - the box can be removed from the holder and quickly brought to the injured person - easy access to the bottles due to spacious design - easy to open - mounted on the wall or standing on the floor - solid/durable material

Additional Information

No.	FG-0138
Manufacturer (Brand)	B-SAFETY
Transport temperature	Room temperature

