



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY® Plum QuickStop 5152 Wound Dressing Set

€11.44
plus VAT &
Shipping

Product Images



Description

QuickStop is a sterile wound dressing set with integrated compresses. Elastic, non-adhesive and skin-friendly. QuickStop is universally applicable as a pressure, cover and support dressing. For heavier bleeding, the smaller wound dressings can also be used as an additional compress. - Sterile wound dressing set with 1 compress 17 x 17 cm and 2 compresses 12 x 7.5 cm - Elastic, non-adhesive and skin-friendly - Packing unit: 6 x 3 pieces - Dimensions (H x W): 114 x 60 mm

Additional Information

No.	FG-0150
Manufacturer (Brand)	B-SAFETY
Transport temperature	Room temperature

