



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY® QuickFix Plaster Dispenser 5500 without filling

€16.37
plus VAT &
Shipping

Product Images



Description

QuickFix plaster dispenser system empty, for filling with plum plaster strips. Incl. fastening material and key for theft protection. With our QuickFix plaster dispenser system, you can fill it yourself and thus adapt it to your specific needs. You can combine all plaster refills from the QuickFix plaster assortment individually according to your wishes. The QuickFix plaster dispenser can be used wherever there is a risk of minor injuries and cuts. Due to its small dimensions, the dispenser can be easily attached to almost any workplace and thus exactly where the need is greatest. The QuickFix can be either screwed on or attached with adhesive strips. Product Features - Dimensions (H x W x D): 133,5 x 232,5 x 33 mm - high hygienic standard due to completely and individually packed plaster strips - optimal protection of the plasters against external contamination - offers immediate and visible access to plasters - for wall mounting directly at the endangered workplace - each plaster is hygienically individually packed and unpacks by itself when pulled out already on one side - the plaster can be applied to the wound with one hand - clear indication when the dispenser is empty - can be refilled quickly - the plasters are CE-marked - built-in shelf for storing small utensils - can also be optionally locked with enclosed key

Additional Information

No.	FG-0157
Manufacturer (Brand)	B-SAFETY
Transport temperature	Room temperature

