




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY® QuickFix plaster dispenser 5505 Alu

€35.25
plus VAT &
Shipping

Product Images



Description

QuickFix Plaster Dispenser 5505 fully equipped with two refill packs of 45 plaster strips each with micronized aluminum wound dressing – improved non-stick effect, fastening material and key for theft protection. For minor injuries, the plaster strips can be quickly and easily removed from the dispenser and are already opened on one side, so that the plaster is only applied to the wound! The dispenser can also be refilled with other types of plasters from the QuickFix range. Plaster with micronized aluminum wound dressing. The non-stick layer of aluminum and PE has a higher non-stick performance than a standard PE layer. The plaster is permeable to air and comfortable on the skin. Product features – PZN: 15570909 – Dimensions (H x W x D): 133.5 x 232.5 x 33 mm – high standard of hygiene due to completely and individually packed plaster strips – optimal protection of the plaster against external contamination – improved non-stick effect – air-permeable – skin-friendly – suitable for sensitive skin

Additional Information

No.	FG-0158
Manufacturer (Brand)	B-SAFETY
Transport temperature	Room temperature

