



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



B-SAFETY® QuickFix Plaster Dispenser 5510 Blood Stopper

€31.92
plus VAT &
Shipping

Product Images



Description

QuickFix Plaster Dispenser 5510 completely equipped with two refill packs of 45 Blood Stopper plasters each with hemostatic effect, fastening material and key for theft protection. For minor injuries, the plaster strips can be quickly and easily removed from the dispenser and are already opened on one side, so that the plaster is only applied to the wound! The dispenser can also be refilled with other plaster types from the QuickFix range. Plaster with haemostatic effect. The alginate fibers between the nonwoven fibers of the wound dressing create an artificial clot (gel) inside it, which stops bleeding within one to two minutes and provides optimal moisture conditions. Product features - PZN: 15743899 - Dimensions (H x W x D): 133.5 x 232.5 x 33 mm - high hygienic standard due to completely and individually packed plaster strips - optimal protection of the plaster against external contamination - hemostatic - with alginate fibers permeable to air - improved protective effect - supports the healing process

Additional Information

No.	FG-0160
Manufacturer (Brand)	B-SAFETY
Transport temperature	Room temperature

