



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® capillary one-way stopcocks, ground and polished, without plug, NS 12.5

**€21.51** As low as  
**€19.15**  
plus VAT &  
Shipping

### Product Images



## Description

---

DURAN® capillary one-way stopcock, ground and polished, without plug, bore 2 mm, NS 12.5 Ground and polished DURAN® one-way stopcocks are stopcocks combined with PTFE stopcocks, thus allowing grease-free operation. Special features and advantages: Ground and polished stopcocks without stopcock plugs Grease-free working due to the polished stopcock surface Pitch 1:10 Glass articles for further processing Application: Installation in glass apparatus and equipment

## Additional Information

---

No.	SD-2446
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

