



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Capillary Disposable Stopcocks, ground, without plug, bore 2 mm, NS 12.5

**€14.26** As low as  
**€12.67**  
**plus VAT &  
Shipping**

### Product Images



## Description

---

DURAN® one-way stopcock, ground, without plug, NS 12.5, bore 2 mm DURAN® one-way stopcocks are designed for use as stopcocks. Valve stopcocks are more suitable for fine dosing. Ground DURAN® one-way stopcocks can be combined with – grooved DURAN® stopcocks – threaded stopcocks – hollow DURAN® stopcocks Special features and advantages: Ground stopcocks without stopcocks Glass articles for further processing Application: Coarse dosing of liquids, installation in glass apparatus and equipment

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | SD-2457          |
| Manufacturer (Brand)  | DURAN            |
| Transport temperature | Room temperature |

