



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**DURAN® capillary one-way stopcocks,
complete with PTFE plug, bore 2.5 mm, NS 14.5**

€32.76 As low as
€29.13
**plus VAT &
Shipping**

Product Images



Description

DURAN® capillary stopcock, with PTFE plug, bore 2.5 mm, NS 14.5 DURAN® stopcocks complete with PTFE stopcock plug allow grease-free working. The 3-part stopcock retainer secures the stopcock plug in the stopcock sleeve. Special features and benefits: Polished stopcock sleeve with stopcock plug made of PTFE for grease-free working Securing of the stopcock plug in the stopcock sleeve by means of 3-part stopcock safety device Pitch 1:10 Glassware for further processing Application: Installation in glass apparatus and equipment

Additional Information

No.	SD-2384
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

