



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® patent stopcocks, complete with SBW plug, bore 4 mm, NS 18.8/50

**€54.67** As low  
as  
**€48.75**  
plus VAT &  
Shipping

### Product Images



## Description

---

DURAN® patent cock, with glass plug, bore 4 mm, NS 18.8 DURAN® patent cocks complete with DURAN® patent cock plug are 2-way cocks often used for the construction of Schlenk lines (vacuum inert gas lines, nitrogen distribution rakes). . The tap plug is secured in the tap sleeve by a 3-part tap retainer. For particularly high leakage rate requirements, we recommend the use of PRODURAN® valves. Special features and benefits: 2-way stopcock with threaded stopcock plug The stopcock plug is secured in the stopcock sleeve via a 3-part stopcock retainer. Glass articles for further processing Application: Suitable for shoring in glass apparatus and for vacuum applications, e.g. for the construction of Schlenk lines (vacuum inert gas lines, nitrogen distribution rakes)

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | SD-2410          |
| Manufacturer (Brand)  | DURAN            |
| Transport temperature | Room temperature |

