



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bohlender® BOLA T-connectors-GL GL 45 Ø 32.5 mm

€148.87
plus VAT &
Shipping

Product Images



Description

BOLA T-connectors-GL T-shaped connector made of PTFE, three connections with GL thread. Universally chemically resistant, the medium only comes into contact with PTFE. Dimensions: Thread GL 45 Through-hole mm 32.5 Length without screw fitting

Additional Information

No.	ER-0867
Manufacturer (Brand)	Bohlender
Transport temperature	Room temperature

