



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Holding device for reaction vessels made of stainless steel, DN 150

€956.56
plus VAT &
Shipping

Product Images



Description

Holding device for reaction vessels, DN 150 Reaction vessels can be operated and temperature-controlled in a variety of ways. To avoid moving the reaction apparatus during the experiment, the reaction vessel is ideally suspended from a stand linkage by means of a suspension device. Cooling and heating media can thus be flexibly fed to the reaction vessel from below and removed again. Special features and advantages: For fixed suspension of the lid and the reaction vessel from a stand rod Application: Holder for lids and reaction vessels, suspension

Additional Information

| | |
|-----------------------|------------------|
| No. | SD-2595 |
| Manufacturer (Brand) | DURAN |
| VGKL number | 132907357 |
| Transport temperature | Room temperature |

