



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## **DURAN® plate holder for glass ceramic laboratory protection plates made of chrome nickel steel 155 x 155 mm**

**€216.60  
plus VAT &  
Shipping**

### **Product Images**



## Description

---

Plate holder, for plate 155 x 155 mm In the laboratory, solutions are often heated. This step can take place using various techniques, e.g. in a water or sand bath, with the aid of a heating hood or classically using a Bunsen burner. If the Bunsen burner is used as the heating source, it is recommended to use a stand or square base and an appropriate glass ceramic plate. The vessel to be heated is placed on the glass ceramic plate, the Bunsen burner underneath. Dirt on the glass ceramic plate can be removed with the help of the cleaning scraper. Special features and advantages: Made of temperature-resistant chrome-nickel steel Application: Tripod plate holder for glass ceramic plates, heating , Suitable for mounting on tripods

## Additional Information

---

No.	SD-2602
Manufacturer (Brand)	DURAN
VGKL number	322907857
Transport temperature	Room temperature

