



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® stainless steel quick release fasteners, with 3 retaining segments, suitable for flange DN 200

€111.74
plus VAT &
Shipping

Product Images



Description

Quick-release closure, DN 200 The quick-release closure closes flat-flange vessels and secures the lid. Slipping or uncontrolled opening is not possible. The closure is made of high-quality stainless steel and features three retaining segments. Special features and advantages: With retaining segments Fast closing due to quick-release closure
Application: Quick-release ring for closure, accessories for plan flange vessels

Additional Information

No.	SD-2593
Manufacturer (Brand)	DURAN
VGKL number	132907161
Transport temperature	Room temperature

