



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle angle beaker 1000 ml

€54.40 As low

as

€48.90

**plus VAT &
Shipping**



Description

Scoop cup

Can be mounted on telescopic rods by means of a snap-on connection.

Order No. 7-8611

Angled cup, 1l

Order No. 7-8612

Pendulum cup, 1l

Automatically swings to the horizontal, suitable for narrow shafts.

Accessories: Telescopic rod (Order No. 7-8607 ff.)

Additional Information

No.	KL-0845
Manufacturer (Brand)	Bürkle
old neolab article no.	7-8611
EAN	4058072267674
Transport temperature	Room temperature
diameter	130 mm
Material	Polypropylene (PP)
Material cup	Polypropylene (PP)
Type liquid collector	Scoop cup Angle cup
Type sampler	Liquid collector
Volume cup	1000 ml
for medium	Liquids

