



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Beakers, tall form, with graduation and spout, 50 ml

**€58.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

DURAN® beakers, tall form, with spout, 50 ml Beakers are cylindrical vessels with a straight side wall. Thanks to this shape, beakers are universally applicable laboratory vessels and are usually used as boiling vessels, as collecting vessels and for preparing, mixing and heating solutions. The large opening of the beaker makes it easy to fill. Beakers are available in different designs and shapes: high and low form, with and without spout, with and without scale. Special features and advantages: Precise spout for clean work Easy-to-read scale with line graduation for measuring the contents Burnt-in and therefore permanently durable white printed image Application: heating, diluting solutions, mixing several components, measuring liquid quantities

## Additional Information

---

No.	SD-0219
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112111617
EAN	4058072183653
Manufacturer EAN	4032051002685
Gross weight	0.042
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.03500 kg
diameter	38 mm
Material	Borosilicate glass 3.3
Volume	50 ml
Country of origin	Germany
sterile	No
Region of origin	Rhineland-Palatinate
Packaging volume	0.00209900 cc
Customs tariff number	70172000
Form	high shape
Height	70 mm

