




neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Usbeck beaker 18/8, graduated, 30 ml capacity

€5.60
plus VAT &
Shipping

Product Images



Description

Ideal replacement for glass or plastic beakers when exceptional material resistance is required.

Easy to clean, autoclavable.
Chemically resistant, graduated.

Material: 18/8 stainless steel

Additional Information

No.	UB-0355
Manufacturer (Brand)	Usbeck
old neoLab article no.	1-1178
EAN	4058072401566
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	51 mm
Color	silver
sterile	No
Height	35 mm
Cups and scales properties	with raised graduation with special chemical resistance graduated suitable for foodstuffs
Material	Stainless steel
Color graduation	silver
Division	10 ml
Material cup	Stainless steel
Type cup	Stainless steel mug
Volume cup	30 ml

