



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm beaker low form with spout, (Boro 3.3), 1000 ml

Special Price

€6.30 was

€8.03 As low as

€5.70

**plus VAT &
Shipping**

Product Images



Description

The RASOTHERM® **beakers**, low form with spout are made of especially heat & chemical resistant borosilicate glass 3.3. With graduation.

All beakers meet the highest quality standards and are manufactured according to the specifications of the current ISO standard 3819.

Thus, the beakers are suitable for many purposes in your laboratory.

Additional Information

No.	1-0159
Manufacturer (Brand)	Rasotherm
EAN	4058072080938
autoclavable at 121°C	Yes
diameter	105 mm
Color	transparent (clear)
sterile	No
Form	low form
Height	145 mm
Cups and scales properties	with spout graduated suitable for foodstuffs
Material	Borosilicate glass
Color graduation	white
Material cup	Borosilicate glass
Type cup	Beakers
Volume cup	1000 ml

